

I-STEP **DETECT**

*Instant Drug & Alcohol
Testing Products*

1413 Pittsburgh Road, Valencia, PA 16059
888-822-7120 (412) 384-5060

Abbreviation	Drug Name	Detection Time (Urine)	Detection Time (Saliva)	Cutoff Level
AMP	Amphetamine	2-4 days	1-3 days	1000 ng/ml
BAR	Barbiturates	1-21 days	1-2 days	300 ng/ml
<u>BUP</u>	Buprenorphine	3-5 days	1-2 days	10 ng/ml
BZO	Benzodiazepines	1-14 days	1-2 days	300 ng/ml
COC	Cocaine	1-7 days	1-2 days	150 ng/ml
COT	Nicotine metabolite	2-4 days	2-4 days	200 ng/ml
EDDP	Methadone metabolite	1-7 days	1-2 days	300 ng/ml
ETG	Ethyl Glucuronide (alcohol)	1-3 days	6-24 hrs	500 ng/ml
FEN	Fentanyl	1-4 days	1-3 days	50 or 100 ng/ml
K2	Spice (synthetic THC)	1-3 days	12-24 hrs	25 or 50 ng/ml
<u>mAMP</u>	Methamphetamine	1-4 days	1-3 days	500 ng/ml
MET	Methamphetamine	1-4 days	1-3 days	500 ng/ml
MDMA	MDMA (ecstasy)	1-3 days	1-2 days	500 ng/ml
MTD	Methadone	1-7 days	1-2 days	300 ng/ml
OPI	Opiates	1-3 days	1-2 days	2000 ng/ml
OPI300	Opiates	1-3 days	1-2 days	300 ng/ml
MOP	Opiates	1-3 days	1-2 days	300 ng/ml
MOR	Opiates	1-3 days	1-2 days	300 ng/ml
OXY	Oxycodone	1-3 days	1-2 days	100 ng/ml
<u>PCP</u>	Phencyclidine	2-13 days	1-2 days	25 ng/ml
PPX	Propoxyphene	1-7 days	1-2 days	300 ng/ml
<u>TCA</u>	Tricyclic Antidepressants	1-21 days	1-2 days	1000 ng/ml
THC	Marijuana	3-30+ days	6-12 hrs	50 ng/ml
TRA	Tramadol	1-3 days	1-2 days	100 ng/ml
6-MAM; 6-AM	Heroin metabolite	1-2 days	8 hrs	10 ng/ml